




THE  
**ULTIMATE**  
**PRE-PACKED PREPARATIVE**  
**COLUMN**  
FOR HPLC AND SFC  
**GUARANTEED!**

Axia PREP LC columns offer:

- Increased Performance
- Groundbreaking Lifetimes
- Optimized Loadability
- Increased Reproducibility



**NEW  
PHASES!**

 **phenomenex**<sup>®</sup>  
...breaking with tradition<sup>SM</sup>



# The Axia™ Advantage

Available in over 40 unique achiral and chiral selectivities, Axia advanced preparative column packings and column hardware designs offer several advantages. Unlike traditional column packing methods, the Axia packing method offers increased sorbent bed density for increased performance and eliminates media bed collapse as a source of premature column failure in preparative HPLC/SFC columns.



guarantee

If Axia packed columns do not provide at least an equivalent separation as compared to a competing preparative column of the same particle size, same phase, and dimensions, return the product with comparative data within 45 days for a FULL REFUND. Only applies to 21.2 mm ID columns.

**Axia™ Technology** ..... pp. 4-7

Award winning column packing technology

**Media Selection Chart** ..... pp. 8-9

Quickly find the media for your purification needs

**Media Selectivity Details** ..... pp. 10-27

Achiral Chemistries

Core-Shell Media

**NEW** Kinetex® ..... pp. 10-13

Aeris™ ..... pp. 14-15

Fully Porous Media

Gemini® ..... pp. 16-17

Synergi™ ..... pp. 18-19

Luna® ..... pp. 20-21

**NEW** Luna® Omega ..... pp. 22-23

Jupiter® ..... p. 24

Clarity® ..... p. 25

Chiral Chemistries

**NEW** Lux® ..... pp. 26-29**PhenoLogix** ..... pp. 30-31

Screening, Optimization and Scale-up Services

**SecurityGuard™ PREP** ..... pp. 32-33

PREP HPLC/SFC Column Protection

**Ordering Information** ..... pp. 33-35

“

I find Axia columns to be very robust and durable. I often use the prep column for much longer than predicted with reproducible peaks. This saves us a significant amount of money.

David Wisnoski  
GlaxoSmithKline, USA

”

# Axia™ Technology vs. Traditional "OBD™" Prep Column Packing


## Axia Packing Technology

Axia packed preparative columns involve a single axial compression step unlike conventional packed preparative columns. The ideal column bed density is custom calculated and automated for each specific media and column size. Computer control of the entire process ensures both proper bed density and column uniformity every time.

During the Axia packing process, the packing piston is locked in place, eliminating any decompression and then recompression of the media sorbent, thus maintaining media and column bed integrity. This solves common lifetime and performance problems associated with conventional packing processes for preparative columns.



## Axia Packing Process Involves: Compression → Final Column



U.S. Patent No. 7, 674, 383

## Traditional Slurry Packing

Traditional slurry packing processes, like the Waters® OBD (Optimum Bed Density) column packing approach, involve the column being removed from the column packing station once it is packed.

Several potential problems with this packing method are:

- Variability in column performance due to increased number of manual operations required for assembly
- Potential silica media damage during recompression
- Level of process control is based on traditional slurry packing technology



## Conventional Packing Process Involves: Compression → Decompression → Recompression → Final Column




Diagram from Waters Corporation U.S. Patent No. 7,399,410


Comparative separations may not be representative of all applications.



## Axia™ packed columns produce uniform media bed with intact particles


The highly tuned patented process and hardware eliminates potential decompression ensuring bed stability and optimal packing density.

The media found on the inlet frit of the Axia packed column shows no signs of damage unlike the media found on inlet frit of traditionally packed prep columns.



## Traditional packed preparative columns produce non-uniform media beds with sheared and crushed particles

Decompression and then recompression during packing can damage the media and lead to increased column-to-column variability, flow disturbances, and decreased column lifetimes.



“

**We are using chromatography media from Phenomenex for GPL/GMP purposes, therefore we audited Phenomenex USA as a manufacturer. From the beginning, we were impressed with Phenomenex and the attitude of their employees. Phenomenex is a unique company in many aspects. Their degree of dedication to customer service, to the organization of the QMS system and last but not least the positive atmosphere in the company is impressive. The outcome of the audit was to our fullest satisfaction.**

”

*Major Generic Pharma Company, Europe*



# Axia™ Technology Outperforms Traditional Packing Processes!

Because of the constant pressure placed on the integrated packing piston, Axia packed columns possess the dynamic capability of maintaining a consistent, homogeneous media bed. This results in superior column performance no matter which media selectivity you choose.


To better understand how much Axia technology improves column performance over traditionally slurry packed preparative columns we scaled-up a 5 µm Lux® Cellulose-1 chiral media analytical column and packed the same media into two different

150 x 21.2 mm I.D. columns. One column was packed using Axia technology and the other prep column was packed using the traditional slurry packing process.


The Axia packing technology had a substantial increase in column efficiency resulting in increased resolution over traditionally packed preparative columns. With increased resolution you are able to increase your sample load enabling you to purify more target compound(s) per purification run. This equates to better throughput and economics.

## Warfarin Chiral Purification in Normal Phase Mode


### Analytical




**Column:** Lux 5 µm Cellulose-1  
**Dimensions:** 150 x 4.6 mm  
**Mobile Phase:** Hexane/Ethanol (75:25)  
**Flow Rate:** 1 mL/min  
**Temperature:** Ambient  
**Inj. Volume:** 100 µL



### Standard Packing and Hardware



### Axia Technology and Hardware



Conditions for both PREP columns:

**Media:** Lux 5 µm Cellulose-1  
**Dimensions:** 150 x 21.2 mm  
**Mobile Phase:** Hexane / Ethanol (75:25)

**Flow Rate:** 20 mL/min  
**Temperature:** Ambient  
**Inj. Volume:** 2 mL

Column (mm)	Analytical 150 x 4.6	Standard 150 x 21.2	Axia 150 x 21.2
Mass Loaded (mg)	2	40	40
Resolution*	1.5	2.85	3.72
Plates (N)	117	535	760

**42% Increase in Efficiency**

\* Resolution calculated with peak width at baseline and center retention time due to the overloaded peaks being off-scale

**Tip:**


For more detailed information on this warfarin application, view application at:

[www.phenomenex.com/tn9002](http://www.phenomenex.com/tn9002)


## Unmatched Column Reproducibility

The completely automated Axia™ packing system provides feedback control and infinite tuning of packing density for specific media characteristics such as mechanical strength and porosity. An optimum higher bed density can be consistently reproduced column-to-column.


This directly translates into consistent efficiency and peak asymmetry measurements and decreases the column variability seen in traditionally packed preparative columns.



Average Efficiency (N) with Synergi™  
4µm Hydro-RP 100 x 21.2mm



Average Peak Asymmetry with Gemini®  
5µm C18 50 x 21.2mm



**9% Increase in Packing, Density, and More Uniformly Packed Bed**



**Axia columns provide me with first rate quality and engineering. Reliability, reproducibility, and durability are provided with all Axia columns that I use. I can literally purify 2500 samples per column. The time and cost savings are tremendous.**

*Derrick Miyao*  
—Large Biotech Manufacturer, USA




**We have used Phenomenex Axia prep-HPLC columns for several years and they consistently provide excellent separation and reproducibility for a variety of different compounds.**

*Jeremy R. Wolf*  
ABC Laboratories, USA




# Phase Selection Chart

Quickly find the media for your purification needs.







**Kinetex**

1st Core-Shell Preparative Column Ever!  
Pages 10-13



**Synergi<sup>TM</sup>**

Unique Chemistries for Complex Mixtures  
Pages 18-19



**Luna**

Proven Purification Performance  
Pages 20-21



**Jupiter**

Increase Loadability for Biomolecule Separations  
Page 24



**Luna Omega**

High Efficiency Polar and Non-Polar Purifications  
Pages 22-23



**Clarity**

Purification of Synthetic Oligonucleotides  
Page 25



**Aeris**

Core-Shell Peptide Media  
Pages 14-15



**Lux**

Polysaccharide Supports with Excellent Enantioselectivity  
Pages 26-29



**Gemini<sup>®</sup>**

High pH Separations  
Pages 16-17

# First Core-Shell Preparative HPLC/SFC Column Ever!

Kinetex® Core-Shell Technology produces increased efficiencies over traditional, fully porous columns, yielding remarkable chromatographic resolution, higher peak capacities, and greater sensitivity, so labs can get even more out of their HPLC analyses!

The benefits of Kinetex Core-Shell Technology include:


- Increased efficiencies over traditional fully porous columns
- Seamless scalability from HPLC/UHPLC to Preparative LC
- Kinetex 5 µm provides better performance than traditional fully porous 5 and 3 µm materials



## High Column Efficiency


Combining 5 µm Kinetex core-shell and Axia™ technologies can provide the highest separation efficiency of any pre-packed preparative HPLC column.

Waters® XBridge® 5 µm C18 Prep OBD™



**Conditions for both columns:**  
**Columns:** Kinetex 5 µm XB-C18 Axia Packed  
 Waters XBridge 5 µm C18 Prep OBD  
**Dimensions:** 150 x 21.2 mm (Kinetex)  
 150 x 19 mm (XBridge)  
**Mobile Phase:** Water/Acetonitrile (50:50)  
**Injection Volume:** 10 µL  
**Flow Rate:** 25 mL/min

Kinetex 5 µm XB-C18 Axia Packed



**Temperature:** Ambient  
**Detection:** UV @ 254 nm  
**Sample:** 1. Uracil  
 2. Acetophenone  
 3. Toluene  
 4. Naphthalene

Key: ● Best Suited ○ Very Good


Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Applications				Type of Compounds				Loading	
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	Available Surface Area
Kinetex C18	1.3, 1.7, 2.6, 5	100	200	12	1.5-8.5*	●	●				○	●	●	○	○
Kinetex XB-C18	1.7, 2.6, 3.5, 5	100	200	10	1.5-8.5*	●	●				○	●	●	○	○
Kinetex EVO C18	1.7, 2.6, 5	100	200	11	1-12	●	●				○	○	●	●	○
Kinetex C8	1.7, 2.6, 5	100	200	8	1.5-8.5*	●	●				○	●	●	○	○
Kinetex Phenyl-Hexyl	1.7, 2.6, 5	100	200	11	1.5-8.5*	●	●				○	●	○	●	○
Kinetex Biphenyl	1.7, 2.6, 5	100	200	11	1.5-8.5*	●	●				○	●	○	●	○
Kinetex HILIC	1.7, 2.6, 5	100	200	0	2.0-7.5	●	○				○	●		○	○
Kinetex F5	1.7, 2.6, 5	100	200	9%	1.5-8.5	●	●				●	●	○	○	○

\*pH stability under gradient conditions. pH stability is 1.5-10.0 under isocratic conditions. Comparative separations may not be representative of all applications.


## Excellent Loadability!

With narrower peak widths than fully porous columns across every sample load, Axia™ packed Kinetex 5 μm columns give you the capability of increased sample load and higher throughput for vastly improved purification performance and economics.

### Waters® XBridge® 5 μm C18 Prep OBD™



### Kinetex 5 μm C18 Axia Packed



**Conditions for both columns:**

**Columns:** Kinetex 5 μm C18 Axia Packed  
 XBridge 5 μm C18 Prep OBD  
**Dimensions:** 50 x 21.2 mm (Kinetex)  
 50 x 19 mm (XBridge)  
**Mobile Phase:** A: Water with 0.5% Formic acid  
 B: Acetonitrile with 0.5% Formic acid  
**Gradient:**

Time (min)	% B
0	20
8	50
11	50

**Flow Rate:** 30 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 254 nm  
**Sample:** 200 mg/mL in DMSO  
 1. Doxepin (From1 - 500 mg on-column)  
 2. Amitriptyline (From1 - 500 mg on-column)



**Kinetex Axia Preparative columns are fantastic! I currently use two Kinetex 5 μm C18 150 x 21.2mm columns in parallel for high throughput purifications (<100mg scale), and Kinetex core-shell media delivers significantly improved peak shape and lower back pressure compared to many of the industry. I can also analyze quickly my purified fractions with the same core-shell phase on my analytical UPLC® system.**



*Chris DeVore*  
 Neurocrine Biosciences, USA

**Tip:**


If you would like to see a loading study performed with the combination of Axia Packing, view application at: [www.phenomenex.com/tn1058](http://www.phenomenex.com/tn1058)

## Seamless Scalability from HPLC/UHPLC to PREP


Kinetex® packed with Axia™ technology makes it the first core-shell sorbent commercially available for small-scale preparative applications. Combine this with the fact that the entire Kinetex

core-shell line is fully scalable from 1.3 µm to 5 µm, this means that transferring high performance HPLC/UHPLC methods to preparative HPLC and SFC formats is fast and simple.

### Waters® XBridge® 5 µm C18 150 x 4.6 mm



### Kinetex 5 µm EVO C18 150 x 4.6 mm




Analytical to PREP Scalability

#### Conditions for all columns:

- Columns:** Kinetex 5 µm EVO C18  
XBridge 5 µm C18
- Dimensions:** 150 x 4.6 mm  
150 x 21.2 mm (Kinetex AXIA Packed)
- Mobile Phase:** A: 0.1 % TFA in Water  
B: 0.1 % TFA in Acetonitrile
- Gradient:** 20 % to 70 % B over 10 minutes
- Flow Rate:** 1.5 mL/min  
30 mL/min (Kinetex AXIA)
- Temperature:** Ambient
- Detection:** UV @ 254 nm
- Sample:** Proprietary Pharmaceutical Sample

### Kinetex 5 µm EVO C18 150 x 21.2 mm AXIA



“

**My Axia packed column has a great efficiency for the separation of several classes of natural compounds. Due to its low back pressure and therefore high flow work conditions, time for conditioning the columns is sped up greatly!**

Sylvian Cretton  
-Europe

”

Comparative separations may not be representative of all applications.


**Tip:**

For more information on the power of Kinetex core-shell scalability, view application at: [www.phenomenex.com/tn1135](http://www.phenomenex.com/tn1135)


# A Broad Spectrum of Column Selectivities

Kinetex® core-shell columns are available in a wide range of stationary phases, allowing you to optimize your separation for maximum resolution and loadability across HPLC, UHPLC, and preparative HPLC and SFC applications.


## C18




## XB-C18




## EVO C18




## C8




## Biphenyl



## Phenyl-Hexyl



## Kinetex 5µm F5



### Conditions for all columns:

**Columns:** Kinetex 5 µm C18  
Kinetex 5 µm XB-C18  
Kinetex 5 µm EVO C18  
Kinetex 5 µm C8  
Kinetex 5 µm Biphenyl  
Kinetex 5 µm Phenyl-Hexyl

**Dimensions:** 100 x 4.6 mm

**Mobile Phase:** A: 0.1% TFA in Water  
B: 0.1% TFA in Acetonitrile

**Gradient:** Time (min) % B  
0 5  
20 20  
22.5 20  
25 5

**Flow Rate:** 1.5 mL/min

**Temperature:** 22 °C

**Detection:** UV @ 330 nm

**Sample:** 1. Chlorogenic Acid  
Others: Antioxidants from green coffee

### Tip:

For more information on Chlorogenic Acids from Green Coffee by HPLC, view application at: [www.phenomenex.com/tn1134](http://www.phenomenex.com/tn1134)




Ordering information on page 34



# Increased Performance for Peptide Purifications


Based on core-shell particle technology, Aeris™ PEPTIDE media is designed with small pores (100 Å), an inert XB-C18 surface chemistry, and multiple particle sizes to meet the selectivity, resolution and loading demands of chemists working with synthetic peptides. The benefits of Aeris PEPTIDE columns include:

- Optimized media for peptide purifications
- Multiple particle size options for method development flexibility and peptide impurity analysis
- Seamless scalability from HPLC/UHPLC to preparative HPLC



XB-C18 chemistry best suited for resolving peptides

## Multiple Particle Sizes For Added Versatility



Packing Material	Particle Size (μm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Applications				Type of Compounds				Loading	
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	Available Surface Area
Aeris PEPTIDE	1.7, 2.6, 3.6, 5	100	200	12	1.5-9.0	●	●					◐	●	◐	◐


# Develop, Purify, and Analyze Peptide Fractions with One Media



Aeris PEPTIDE is fully scalable in retention and selectivity with its 4 unique particle sizes (1.7 μm, 2.6 μm, 3.6 μm, and 5 μm) for easy transfer from HPLC and UHPLC methods to preparative applications.

## Aeris PEPTIDE 2.6 μm XB-C18

### Analytical method




**Column:** Aeris PEPTIDE 2.6 μm XB-C18  
**Dimensions:** 150 x 4.6 mm  
**Part No.:** [00F-4505-E0](#)  
**Injection Volume:** 10 μL  
**Flow Rate:** 1 mL/min  
**Sample:** Crude peptide mix



## Aeris PEPTIDE 5 μm XB-C18

### Preparative scale-up and fraction collection




**Column:** Aeris PEPTIDE 5 μm XB-C18  
 Axia Packed  
**Dimensions:** 150 x 21.2 mm  
**Part No.:** [00F-4632-P0-AX](#)  
**Injection Volume:** 1 mL  
**Flow Rate:** 20 mL/min  
**Sample:** Crude peptide mix



## Aeris PEPTIDE 2.6 μm XB-C18

### Analytical fraction analysis



**Column:** Aeris PEPTIDE 2.6 μm XB-C18  
**Dimensions:** 150 x 4.6 mm  
**Part No.:** [00F-4505-E0](#)  
**Injection Volume:** 10 μL  
**Flow Rate:** 1 mL/min  
**Sample:** Purified Fractions



**Conditions for all separations (except as noted):**  
**Mobile Phase:** A: 0.1 % TFA in Water  
 B: 0.1 % TFA in Acetonitrile  
**Gradient:** Linear 85:15 (A/B) to 5:95 (A/B) over 10 minutes  
**Temperature:** Ambient  
**Detection:** UV @ 210 nm



Ordering information on page 34

# Setting the Standard for pH Method Development



Gemini features a pH stability from 1-12, making it optimal for high alkaline washes and high pH purifications of basic drugs.

### Optimized parameters include:


- Innovative surface layer for increased pH stability
- High-surface area for increased loading
- Silica smoothness for stable packing beds
- Bonding density for excellent reproducibility

## Second-Generation TWIN-NX™ Technology

### Gemini NX-C18


TWIN-NX technology uses an improved, patented organo-silica grafting process which incorporates highly stabilizing ethane cross-linking. These organic groups are evenly incorporated into the grafted layers on the silica surface while maintaining a pure silica core. This not only provides resistance to high pH attack, but also maintains the high efficiency and mechanical strength of a silica particle.

\*This bonding technology is also available in Core-Shell media. See Kinetex® EVO on page 12.




Dramatically improve sample resolution, productivity and performance of any preparative column media with Axia™ column hardware and packing technology. Axia packed prep columns offer the opportunity for longer lifetime, higher loading and increased throughput.

### Gemini 5 µm NX-C18 Axia Packed



### Waters® XBridge® 5 µm C18 Prep OBD™



Conditions for both columns:

**Column:** Gemini 5 µm NX-C18  
Waters 5 µm XBridge  
**Dimensions:** 150 x 21.2 mm (Gemini)  
150 x 19 mm (XBridge)  
**Mobile Phase:** A: 20 mM Ammonium bicarbonate pH 10.0  
B: Acetonitrile

**Gradient:** Time (min) % B  
0.1 50  
5.1 95  
6 95  
6.5 50  
8.9 50

**Flow Rate:** 25 mL/min  
**Temperature:** 22 °C  
**Detection:** UV @ 268 nm  
**Sample:** 1. Reserpine  
2. Unknown

Key: ● Best Suited ○ Very Good


Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Applications					Type of Compounds				Loading	
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	Available Surface Area	
Gemini C18	3, 5, 10	110	375	14	1-12	●	●				●	○	●	●	●	●
Gemini C6-Phenyl	3, 5	110	375	12	1-12	●	●				○		●	○	●	●
Gemini NX-C18	3, 5, 10	110	375	14	1-12	●	●				●		●	●	●	●

Comparative separations may not be representative of all applications.

# Flexibility in pH Adjustments Allows for Increased Purification Performance

Separating basic compounds at higher pH levels produces dramatic changes when compared to low pH conditions. At pH 10.5, the basic compounds become neutralized and are more hydrophobic. The retention for the uncharged basic compounds increases providing an increase in separation along with superior peak shapes.


## Gemini® NX-C18 with 0.5% TFA



**Column:** Gemini NX-C18 5 µm  
**Dimensions:** 50 x 21.2 mm  
**Mobile Phase:** A: 0.5% TFA in Water  
 B: Acetonitrile  
**Gradient:** 5% B to 95% B in 5 min  
**Flow Rate:** 30 mL/min  
**Detection:** UV @ 254 nm  
**Sample:** 1. Diphenhydramine  
 2. Propranolol


Excellent resolution at high pH

## Gemini NX-C18 with 0.2% NH<sub>4</sub>OH



**Column:** Gemini NX-C18 5 µm  
**Dimensions:** 50 x 21.2 mm  
**Mobile Phase:** A: 0.2% NH<sub>4</sub>OH in Water  
 B: Acetonitrile  
**Gradient:** 5% B to 95% B in 5 min  
**Flow Rate:** 30 mL/min  
**Detection:** UV @ 254 nm  
**Sample:** 1. Diphenhydramine  
 2. Propranolol


## Gemini NX-C18 with 0.5% TFA



**Column:** Gemini NX-C18 5 µm  
**Dimensions:** 50 x 21.2 mm  
**Mobile Phase:** A: 0.5% TFA in Water  
 B: Acetonitrile  
**Gradient:** 5% B to 95% B in 5 min  
**Flow Rate:** 30 mL/min  
**Detection:** UV @ 254 nm  
**Sample:** 1. Diphenhydramine  
 2. Oxybutynin  
 3. Terfenadine

Separation shape improvement provides opportunity for increased loading

## Gemini NX-C18 with 0.2% NH<sub>4</sub>OH



**Column:** Gemini NX-C18 5 µm  
**Dimensions:** 50 x 21.2 mm  
**Mobile Phase:** A: 0.2% NH<sub>4</sub>OH in Water  
 B: Acetonitrile  
**Gradient:** 5% B to 95% B in 5 min  
**Flow Rate:** 30 mL/min  
**Detection:** UV @ 254 nm  
**Sample:** 1. Diphenhydramine  
 2. Oxybutynin  
 3. Terfenadine



Our Phenomenex Gemini and Luna® Axia™ packed columns are the work-horses in our lab. These columns exhibit outstanding performance for challenging separations while also handling a high workload for standard separations. Longevity has also been excellent with some columns lasting 2 years or more. Dependability is so important in my line of work and these columns never disappoint!!



-Major Pharmaceutical Company, USA



Ordering information on page 34

### Tip:

If you want longer Gemini NX-C18 Axia packed column lifetimes, view a lifetime study application at: [www.phenomenex.com/tn1138](http://www.phenomenex.com/tn1138)


# Increased Loading with Unique Selectivities

Synergi is available in four unique phases, each offering dramatic differences in:


- Selectivity
- Retention time
- Resolution

The unique selectivity profiles found within the Synergi product line offer complementary selectivity to the standard C18, C8, or silica phases traditionally employed in preparative HPLC.


## Synergi Polar-RP For Polar and Aromatic Mixtures (100% Aqueous Stable)




## Synergi Fusion-RP Balanced Non-polar and Polar Performance (100% Aqueous Stable)



## Synergi Hydro-RP Strong Non-polar and Polar Retention (100% Aqueous Stable)



## Synergi Max-RP Excellent for Basic Compounds at Neutral pH



Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Applications					Type of Compounds				Loading Available Surface Area
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	
Synergi Fusion-RP	4, 10	80	475	12	1.5-9.0*	●	●				●	●	●	●	●
Synergi Max-RP	4, 10	80	475	17	1.5-9.0*	●	●				●		●	●	●
Synergi Hydro-RP	4, 10	80	475	19	1.5-7.5	●	●				●	●	●	●	●
Synergi Polar-RP	4, 10	80	475	11	1.5-7	●	●				●	●		●	●



\*pH stability under gradient conditions. pH stability is 1.5-10.0 under isocratic conditions.



## Selectivity Like No Other

Offering a balanced combination of hydrophobic and polar selectivity, Synergi™ Fusion-RP separates compounds exhibiting moderately polar and hydrophobic characteristics.

### Hydrophobic basic compounds







**Conditions for all columns:**

- Columns:** Synergi 4 µm Fusion-RP  
Typical C18
- Dimensions:** 150 x 4.6 mm
- Mobile Phase:** 20 mM Potassium Phosphate, pH 2.5 / Acetonitrile (75:25)
- Flow Rate:** 1.0 mL/min
- Detection:** UV @ 210 nm
- Sample:** 1. Maleic acid  
2. Chlorpheniramine  
3. Triprolidine  
4. Diphenhydramine

The slightest variations in compound polarity and aromaticity are exploited by Synergi Polar-RP to achieve the greatest separation between polar and/or aromatic compounds.

### Increased resolution of polar compounds with Synergi Polar-RP compared to traditional C18 phases




**Conditions for all columns:**

- Columns:** Synergi 4 µm Polar-RP  
Waters 5 µm SymmetryShield™ RPC18  
Waters 5 µm Symmetry® C18  
Waters 5 µm XTerra® RP18
- Dimensions:** 150 x 4.6 mm
- Mobile Phase:** 20 mM Potassium phosphate pH 3 / Methanol (50:50)
- Flow Rate:** 1.0 mL/min
- Detection:** UV @ 230 nm
- Temperature:** Ambient
- Injection:** 2 µL
- Sample:** 1. Metaproterenol (0.4 µg)      4. Alprenolol (0.3 µg)  
2. Pindolol (0.6 µg)                      5. Propranolol (0.04 µg)  
3. Metoprolol (0.15 µg)                  6. Ethylparaben (0.4 µg)

“ We regularly use RP stationary phases from Phenomenex for our separation problems. Especially Synergi Polar-RP which was found to often show the desired selectivity, distinguishing this phase from other RP phases. ”

*CARBOGEN AMCIS, Switzerland*

Comparative separations may not be representative of all applications.




# Media for One of the World's Leading PREP HPLC Columns

Luna® high surface area (400 m<sup>2</sup>/g) silica packing materials provide optimized parameters specifically designed for the purification of small molecules and peptides. This media allows high loading with excellent lifetimes.

## Optimized loading parameters include:

- Silica smoothness for stable packed beds
- Optimum pore size/distribution provide outstanding performance
- High pore volume offers increased surface area
- Fine tuned bonding density for excellent reproducibility
- Greater loading capacity with an extended pH range of 1.5 to 10.0\*



We routinely use Axia™ packed columns from Phenomenex for peptide purifications. Among various preparative HPLC columns we have used, the Axia packed Luna columns (5 µm) stand out. We have been very satisfied with the increased loading capacity and excellent performance.



Guangcheng Jiang  
Ferring Research Institute, Inc., USA

Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Applications					Type of Compounds				Loading
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	
Luna C18(2)	3, 5, 10, 10-PREP, 15	100	400	17.5	1.5-9.0*	●	●				●		●	●	●
Luna C8(2)	3, 5, 10, 10-PREP, 15	100	400	13.5	1.5-9.0*	●	●				●		●	●	●
Luna C5	5, 10	100	440	12.5	1.5-9.0*	●	●				●		●	●	●
Luna Phenyl-Hexyl	3, 5, 10, 10-PREP, 15	100	400	17.5	1.5-9.0*	●	●				●	●	●	●	●
Luna Silica(2)	3, 5, 10, 10-PREP, 15	100	400	-	-	●					●	●		●	●
Luna CN	3, 5, 10	100	400	7	1.5-7.0	●					●	●		●	●
Luna NH <sub>2</sub>	3, 5, 10	100	400	9.5	1.5-11.0	●					●	●		●	●
Luna SCX	5, 10	100	400		2-7	●						●		●	●
Luna HILIC	3, 5	200	200	5.7	1.5-8.0	●	●				●	●		●	●
Luna PFP(2)	5, 10	100	400	11.5	1.5-8.0	●						●		●	●


\*pH stability under gradient conditions. pH stability is 1.5-10.0 under isocratic conditions.



## Simple Scale-Up

Axia™ column technology provides the same high efficiency chromatographic performance for preparative scale columns (21.2, 30, and 50 mm ID) as obtained in 4.6 mm ID analytical columns. This improvement in preparative column performance across

all lengths and internal diameters makes it easier to select the appropriate column size to achieve the desired purity and yield without having to compromise on performance.




## Proven Media for Peptide Purifications

An optimal compromise between throughput, recovery, and yield. Perform high loading (0.74g on column) and achieve high purity (>98 %) in a single purification run.

### Preparative Purification of Bivalirudin


(20 amino acid peptide also known as Angiomax®)

#### Purification Elution Profile at 1.5 % Specific Load



**Column:** Luna 10 µm-PREP C8(2)  
**Dimensions:** 250 x 21.2 mm  
**Part No.:** 00G-4323-PO-AX  
**Mobile Phase:** A: 100 mM Ammonium acetate pH 4.7 in Water  
B: Acetonitrile  
**Gradient:** 10 to 50 % B in 40 min; hold at 80 % B for 5 min; re-equilibration at 10 % B for 10 min  
**Flow Rate:** 21 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 280 nm  
**Injection Volume:** 105 mL  
**Sample Concentration:** 7 mg/mL in water  
**Sample:** Crude Bivalirudin in Water

#### Purity Confirmation of Combined Fractions



**11 Combined fractions 27.8 – 29.8 min;**  
**Recovery 80.5 % with purity ≥ 98.5 %**

Peak No.	Time (min)	Area	Area %
1	12.74	73.7	0.35
2	13.83	40.6	0.19
3	14.37	21118.7	98.53
4	15.858	200.5	0.93

**Column:** Luna 5 µm C8(2)  
**Dimensions:** 250 x 4.6 mm  
**Part No.:** 00G-4249-E0  
**Mobile Phase:** A: 0.1 % TFA in Water  
B: 0.1 % TFA in Acetonitrile  
**Gradient:** 20 % to 50 % B in 30 min

**Flow Rate:** 1 mL/min  
**Temperature:** 25 °C  
**Detection:** UV @ 220 nm  
**Injection Volume:** 2 µL  
**Sample:** Combined Fractions



Ordering information on page 35


# Luna<sup>®</sup> Omega **NEW**

## Cutting Edge Fully Porous Silica Particle

Luna is one of the most recognized HPLC brands on the market, delivering high efficiency, ruggedness, reproducibility, and dependability for a wide range of analyses. The new Luna Omega builds upon this legacy with an innovative yet rugged silica particle architecture, designed and manufactured by Phenomenex based on more than 20 years of applied knowledge, invention, and customer experience.

## Novel Design and Manufacturing Process

Within the novel manufacturing process of Luna Omega silica, we implement a proprietary processing technique to gain greater particle inertness, a stronger particle morphology, and more consistent porosity.




Packing Material	Particle Size (μm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Applications					Type of Compounds				Loading
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	
<b>NEW</b> Luna Omega Polar C18	1.6, 3, 5	100	260	9%	1.5-8.5*	●	●				●	●	●	●	●
<b>NEW</b> Luna Omega PS C18	1.6, 3, 5	100	260	9%	1.5-8.5*	●	●				●	●	●	●	●

\*pH stability under gradient conditions. pH stability is 1.5-10 under isocratic conditions.


## UHPLC to HPLC to PREP Scalability

With direct selectivity scalability from Luna Omega 1.6 µm to 5 µm you can fluidly transfer methods from UHPLC platforms to HPLC and preparative instrumentation. Additionally, you can easily go in reverse and use a Luna Omega 1.6 µm to analyze fractions taken from a Luna Omega 5 µm preparative column.

### Luna Omega Phase Selection




100 % aqueous stability and enhanced selectivity/retention for polar analytes without diminishing useful non-polar retention. The C18 ligand provides general hydrophobic interactions while a polar modified particle surface provides enhanced polar compound retention.



Unique, 100 % aqueous stable mixed-mode phase that provides both polar and non-polar retention. The surface contains a positive charged ligand which aids in the retention of acidic compounds through ionic interactions, while the C18 ligand promotes general reversed phase retention. The positively charged surface also improves basic compound peaks shape through ionic repulsion.

### Direct Scalability 1.6 µm to 5 µm




Conditions for all columns:

- Columns:** Luna Omega 1.6 µm Polar C18  
Luna Omega 5 µm Polar C18
- Dimension:** 50 x 2.1 mm
- Mobile Phase:** A: Water with 0.1 % Formic Acid  
B: Acetonitrile with 0.1 % Formic Acid
- Gradient:**

Time (min)	% B
0	5
5	95
- Flow Rate:** 0.4 mL/min
- Temperature:** 30 °C
- Detection:** UV @ 254 nm
- Sample:**
  1. Uracil
  2. Pindolol
  3. Chlorpheniramine
  4. Nortriptyline
  5. 3-Methyl-4-nitrobenzoic acid
  6. 5-Methyl salicylaldehyde
  7. Hexanophenone


### UHPLC to HPLC to PREP

#### UHPLC Scout




Dimension: 50 x 2.1 mm  
Flow Rate: 0.2 mL/min

#### Optimized HPLC



Dimension: 150 x 4.6 mm  
Flow Rate: 1.5 mL/min

#### Preparative Purification



Dimension: 150 x 21.2 mm  
Flow Rate: 30 mL/min

Conditions for all columns (as noted):

- Columns:** Luna Omega 5 µm PS C18
- Mobile Phase:** A: Water with 0.1 % TFA  
B: Acetonitrile with 0.1 % TFA
- Gradient:**

Time (min)	% B
0	10
15	90
- Temperature:** 22 °C
- Detection:** UV @ 254 nm
- Sample:**
  1. Impurity
  2. Proprietary API
  3. Impurity

Analyze Prep Fractions  
via UHPLC



Ordering information on page 35



# Media for Biomolecules







The Jupiter HPLC column portfolio, including Jupiter 300 and Jupiter Proteo, offers optimized reversed phase solutions for peptide and protein purification. Identify, purify, and analyze almost any protein with Jupiter columns.

## Jupiter Proteo 90 Å

- For separation of proteins and peptides < 10,000 MW
- C12 bonded onto an ultra-high surface area (475 m<sup>2</sup>/g) silica for increased peak capacity and resolution of peptide separations
- Direct scale-up from analytical to preparative and bulk materials

### Resolve Peptides with Similar Hydrophobicity

Jupiter Proteo is able to fully resolve peptides that differ in hydrophobicity by one methyl group.



**Columns:** Phenomenex Jupiter 4 µm Proteo 90 Å  
 Agilent Technologies ZORBAX 5 µm SB-C18 300 Å  
 Alltech Associates, Inc. Vydac 5 µm MS54 300 Å  
 Alltech Associates, Inc. Vydac 5 µm TP54 300 Å

**Dimensions:** 250 x 4.6 mm

**Mobile Phase:** A: 0.1% TFA in Water  
 B: 0.085% TFA in Acetonitrile

**Gradient:** A/B (95:5) to A/B (55:45) in 20 minutes

**Flow Rate:** 1 mL/min

**Temperature:** 22 °C

**Detection:** UV @ 214 nm


**Sample:** 1. NH<sub>2</sub>-Arg-Gly-Gly-Ala-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide  
 2. Ac-Arg-Gly-Gly-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide  
 3. Ac-Arg-Gly-Ala-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide  
 4. Ac-Arg-Gly-Val-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide  
 5. Ac-Arg-Gly-Val-Val-Gly-Gly-Gly-Leu-Gly-Leu-Gly-Lys-Amide

## Jupiter 300 Å

- For separation of proteins > 10,000 MW
- Available with C18 and C4 bonded phases
- 1.5 – 10 pH stability for method ruggedness and easy protein removal
- Direct scale up to preparative and bulk materials

### Compare PEGylated vs. Native Forms of Proteins

Reversed phase separation of PEGylated and native proteins on a Jupiter 300 C4 column. Note the good resolution of multiple PEGylated forms for all proteins tested.



**Columns:** Jupiter 300 5 µm C4 300 Å  
**Dimensions:** 150 x 4.6 mm  
**Part No.:** 00F-4167-E0  
**Mobile Phase:** A: 2% Acetonitrile / 0.1% TFA in Water  
 B: 70% Acetonitrile / 20% IPA / 0.08% TFA in Water  
**Gradient:** A/B (85:15) to A/B (30:70) in 25 min  
**Flow Rate:** 1 mL/min  
**Temperature:** 45 °C  
**Detection:** UV @ 214 nm  
**Sample:** PEGylated and Native Proteins

“ We purchased the Jupiter 300 C18 300 Å column a few months ago and have been quite impressed with its performance. The Jupiter 300 column provides better separation of the proteins. As for reproducibility, the control profiles have not changed since day one of its use. ”

*Major Biotech Company, Europe*

Key: ● Best Suited ○ Very Good

Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Applications					Type of Compounds				Loading
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	
Jupiter C18	5, 10, 15	300	170	13.3	1.5-10.0			●				○	●	○	●
Jupiter C4	5, 10, 15	300	170	5	1.5-10.0			●				○	●	○	●
Jupiter Proteo	4, 10	90	475	15	1.5-10.0	●	●					○	●	○	●

Comparative separations may not be representative of all applications.

# Purification of Synthetic Oligonucleotides




## Clarity Oligo-RP™

Unique media specifically designed for reversed phase purification of oligonucleotides with balanced hydrophobicity and polar selectivity. The media is based on composite particle TWIN™ technology that provides improved selectivity and efficiency for oligonucleotides when compared to competing hybrid, polymer, and silica media.

### RP-HPLC Preparative Purification

- Easily separate N-1 failure sequences from target oligo with >90 % purities
- Purify oligos up to 60 nt in length
- Trityl-off purification of DNA, RNA, thioates, and modified/ labeled oligonucleotides
- 3 µm, 5 µm, 10 µm particles for seamless scaling

### Preparative 20nt DNA Oligo-RP Purification



**Column:** Clarity 3 µm Oligo-RP C18  
**Dimensions:** 50 x 10 mm  
**Part No.:** OOB-4441-NO  
**Mobile Phase:** A: 50 mM TEAA pH 7.5/ 5 % Acetonitrile  
 B: Methanol  
**Gradient:** 10 % to 60 % B in 20 minutes  
**Flow Rate:** 5 mL /min  
**Detection:** UV @ 260 nm  
**Sample:** 20nt DNA

## Clarity Oligo-XT **NEW**

Clarity Oligo-XT, C18 columns have been designed to provide rugged high performance for the LC/MS characterization of synthetic DNA and RNA samples, alongside purification of these targeted oligos. With high efficiency levels from the novel core-shell particle design, this new media provides the necessary separation power to accurately resolve closely related oligo sequences.

- Novel core-shell particle technology with rugged pH stability from 1-12
- 5 µm particles provide extremely low pressure HPLC and Preparative purification solutions
- Seamless scalability between all three particle sizes (1.7 µm, 2.6 µm and 5 µm)

## Clarity Oligo-WAX™

Clarity Oligo-WAX is a crosslinked weak anion-exchange media designed for successful ion-exchange purification of synthetic DNA/RNA. Oligo-WAX is an advantageous combination of purity, capacity, mechanical strength, cost, and efficiency.

- Excellent efficiency column results in > 90 % purities due to good fractionation of closely eluting compounds
- High loading capacity due to very high density ligand
- Increase productivity by running at higher flow rates and pressures



**We have used the Axia™ prep columns and have not had problems with them. I have never had to adjust for retention gaps. This speaks directly to the quality of Phenomenex's phases and the quality of their PREP columns.**



-Major Biotech Company, USA

Packing Material	Particle Size (µm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Applications				Type of Compounds				Loading	
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic		Bases
Oligo-RP	3, 5, 10	110	375	14	1-12	●	●	●	●	●	●	●	●	●	●
Oligo-Wax	10	360	150	-	1-11	●	●	●	●	●	●	●	●	●	●
<b>NEW</b> Oligo-XT	1.7, 2.6, 5	100	200	11	1-12	●	●	●	●	●	●	●	●	●	●

Key: ● Best Suited ○ Very Good

# Complete Chiral Solutions




Achieving optimal chiral separation is easier than ever with seven unique Lux polysaccharide stationary phases to screen. Choose a phase, then transfer the method to lab scale, process, pilot, and commercial scale.

## Lux chiral preparative columns simplify the separation process:

- Unique and traditional phases that increase the success rate of the chiral screen
- Consistent particle size distribution so performance is maintained
- Mechanically strong media for increased stability
- Available in multiple particle sizes for direct scale up (3 μm and 5 μm packed columns for screening and small scale purifications; 10 μm and 20 μm bulk media for process scale purifications)


## Resolve Your Enantiomers with Seven Unique Phases

The Lux family of bulk cellulose and amylose chiral selectors provides a variety of complementary selectivities.




### Lux Cellulose-1

Cellulose tris(3,5-dimethylphenylcarbamate)  
**Guaranteed Alternative to** CHIRALCEL® OD®, OD-H®, OD-3, OD-RH®, and OD-3R




### Lux Cellulose-3

Cellulose tris(4-methylbenzoate)  
**Guaranteed Alternative to** CHIRALCEL OJ®, OJ-H®, OJ-3, OJ-RH®, and OJ-3R




### Lux Amylose-1

Amylose tris(3,5-dimethylphenylcarbamate)  
**Guaranteed Alternative to** CHIRALPAK AD®, AD-H®, AD-3, AD-RH®, and AD-3R




### Lux i-Cellulose-5

Cellulose tris(3,5-dichlorophenylcarbamate)  
**Guaranteed Alternative to** CHIRALPAK® IC and IC-3




### Lux Cellulose-2

Cellulose tris(3-chloro-4-methylphenylcarbamate)  
**Guaranteed Alternative to** CHIRALCEL OZ, OZ-H®, OZ-3, OZ-RH, and OZ-3R



### Lux Cellulose-4

Cellulose tris(4-chloro-3-methylphenylcarbamate)  
**Guaranteed Alternative to** CHIRALCEL OX-H, OX-3, OX-RH, and OX-3R



### Lux Amylose-2

Amylose tris(5-chloro-2-methylphenylcarbamate)  
**Guaranteed Alternative to** CHIRALPAK AY®, AY-H®, AY-3, AY-RH, and AY-3R

Key: ● Best Suited ○ Very Good

Packing Material	Particle Size (μm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Range	Chiral Applications				Type of Chiral Compounds				Loading		
						Small Molecules	Peptides	Proteins	Chiral	Oligonucleotides	Acids	Polar	Hydrophobic	Bases	Available Surface Area	
<b>NEW</b> Lux i-Cellulose-5	3, 5	1,000	-	-	2-9	●			●		○	○	○	○	○	●
Lux Cellulose-1	3, 5, 10, 20	1,000	-	-	2-9	●			●		○	○	○	○	○	●
Lux Cellulose-2	3, 5, 10, 20	1,000	-	-	2-9	●			●		○	○	○	○	○	●
Lux Cellulose-3	3, 5, 10, 20	1,000	-	-	2-9	●			●		○	○	○	○	○	●
Lux Cellulose-4	3, 5, 10, 20	1,000	-	-	2-9	●			●		○	○	○	○	○	●
Lux Amylose-1	3, 5	1,000	-	-	2-9	●			●		○	○	○	○	○	●
Lux Amylose-2	3, 5	1,000	-	-	2-9	●			●		○	○	○	○	○	●

# Column Screening for Optimal Chiral Resolution


Being able to utilize differences in selectivity in each of the seven Lux® columns can help develop methods more efficiently by offering broad and complementary chiral recognition abilities.

In the example below, a simple screen determined which column gave the best separation.


## Etozolin

Based on a five phase screen under reversed phase conditions, the optimal chiral stationary phase for resolving Etozolin is Lux Cellulose-3.


### Lux 5µm Cellulose-3




**Conditions for all columns:**  
**Column:** As noted  
**Dimension:** 250 x 4.6 mm  
**Mobile Phase:** Acetonitrile /  
 20 mM Ammonium  
 bicarbonate with 0.1 % Diethylamine (60:40)  
**Flow Rate:** 1 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 220 nm




### Lux 5µm Amylose-2




### Lux 5µm Cellulose-2



### Lux 5µm Cellulose-1




### Lux 3µm Cellulose-4




## Innovative chiral selector will succeed where others fail


### Lux 5µm Cellulose-4




**Conditions for all columns:**  
**Dimensions:** 250 x 4.6 mm  
**Mobile Phase:** 0.1 % Diethylamine in Hexane /  
 0.1 % Diethylamine in Isopropanol (90:10)  
**Flow Rate:** 1 mL/min  
**Detection:** UV @ 220 nm  
**Temperature:** Ambient



### CHIRALCEL® 5µm OD-H®



### CHIRALPAK® 5µm AD-H®




Columns used for comparison were manufactured by DAICEL Corporation.  
 Phenomenex is in no way affiliated with DAICEL Corporation. Comparative separations may not be representative of all applications.


## Load More with an Increase in Column Length

Axia™ column technology allows separations to scale up directly based on column length. With the 100mm length column a 32 mg/load separation was achieved and an increased sample


load of 80 mg/load was achieved on the longer 250 mm length column. As expected when increasing the load, the peak width and tailing increased but there was no loss of resolution.



No resolution loss with increased sample load




2.5x Load Increase



### Conditions for all columns:

**Columns:** Lux® 5  $\mu$ m Cellulose-1  
**Dimensions:** as noted  
**Mobile Phase:** Methanol / Isopropanol (90:10)  
**Flow Rate:** as noted  
**Detection:** as noted  
**Sample:** Dissolved in mobile phase as noted

**Dimensions:** 100 x 4.6 mm  
**Flow Rate:** 1 mL/min  
**Detection:** UV @ 220 nm  
**Sample:** 5  $\mu$ g in 2  $\mu$ L



Methocarbamol

**Dimensions:** 100 x 21.2 mm  
**Flow Rate:** 20 mL/min  
**Detection:** UV @ 220 nm and 254 nm  
**Sample:** 32 mg in 640  $\mu$ L

**Dimensions:** 250 x 21.2 mm  
**Flow Rate:** 20 mL/min  
**Detection:** UV @ 220 nm and 254 nm  
**Sample:** 80 mg in 1600  $\mu$ L

“

Lux Axia preparative columns are wonderful! I regularly use Lux chiral stationary phases Cellulose-2 and Cellulose-4 and less frequently, the Lux Amylose-2. In our community of chiral analysis/purification scientists, there are some who use the CC4 column instead of the \*equivalent\* Lux Cellulose-4. On several occasions we've seen separation and good peak shape on the Lux Cellulose-4 that was completely missing from the CC4. Customer support and delivery times are always within a few days.

Julia G. Christie  
 GlaxoSmithKline, USA


”



# Easy SFC Scale-up

## SFC Purification of Terfenadine

### Baseline Separation of Enantiomers




Lux® Cellulose-1 offers great peak shape at 220 nm


**Dimensions:** 250 x 4.6 mm  
**Flow Rate:** 2.5 mL/min  
**Detection:** UV @ 220 nm  
**Load:** 300 µg in 10 µL



Analytical and Axia™ packed columns have been extensively tested on various SFC systems and all column ID's and lengths are SFC compatible.



5x Load Increase



Overloading study with increased analytical load showing impurities eluting after major enantiomers only detected at 254 nm

**Dimensions:** 250 x 4.6 mm  
**Flow Rate:** 2.5 mL/min  
**Detection:** UV @ 254 nm  
**Load:** 1.5 mg in 50 µL

**Conditions for all columns:**

**Columns:** Lux 5 µm Cellulose-1

**Mobile Phase:** Methanol with 0.1 % DEA/  
Carbon Dioxide (25:75)

**Column Temperature:** 35 °C

**Polarimeter:** ALP-PDR-Chiral

**Sample:** Terfenadine with ethanol  
dissolution solvent

70x Load Increase



7.5 cycles/hr  
787 mg/hr

High loading capacity media along with stacking injections allow for increased yields and productivity

Closer stacked injections can not be used due to the impurities eluting after the major enantiomers

**Dimensions:** 250 x 21.2 mm  
**Flow Rate:** 50 mL/min  
**Detection:** UV @ 220 nm  
**Load:** 105 mg in 3.5 mL

Tip:

For SFC column screening, use Lux 150 x 3.0 mm ID columns.



Ordering information on page 35

# A New Era of Technical Support Services

## Let Us Do the Work for You

PhenoLogix, our in-house application support lab, saves you time and money by screening multiple scout columns and solvent strategies for new purification methods or revalidating your current methods. We work together to make you successful by minimizing your process purification development time and optimizing your purification method.

### Chiral Screening

- Normal Phase
- Reversed Phase
- Polar Organic
- SFC

### Method Optimization Services

- Fast Turnaround
- Easy Method Transfer
- Continued Support

### Preparative and Process Scale-Up

- Media Screening
- Small Scale Purification
- DAC Packing Assistance




## Your Method, Our Scientists

Quality Products, Advanced Performance, Complete Support

For more information or to begin a project today, please contact your local Phenomenex representative or email us at [phenologix@phenomenex.com](mailto:phenologix@phenomenex.com)

You can also visit us online:

[www.phenomenex.com/phenologix](http://www.phenomenex.com/phenologix)



“

**Our scientists at American Peptide have taken advantage of Phenomenex's column packing services, application development, and project-specific consultation services for some of our most challenging separations.**

”

*American Peptide Company, USA*



# PREP HPLC/SFC Column Protection

## SecurityGuard™ PREP

- Extends preparative column lifetime by as much as 5 times
- Protects columns from samples that precipitate out of solution
- Protects columns from contaminants
- Stable and leak-free up to 60 mL/min within specified pressure ratings


The SecurityGuard PREP system was designed to effectively and inexpensively protect your valuable prep columns from the damaging effects of mobile phase and sample chemical contaminants and particulates, without altering your chromatographic results.

## Lower Your Cost Per Injection!

SecurityGuard isn't only about column protection, it's about lowering your cost per injection! When you increase the number of injections from a single preparative column you're lowering your overall cost per injection. With SecurityGuard PREP, the inexpensive cartridge was ruined while the integrity of the prep column was maintained and its performance restored.


### Forced Degradation Study

**Injection 1:** Axia™ Packed column with SecurityGuard PREP cartridge



**Column:** Luna® 10µm C18(2) Axia Packed  
**Dimension:** 50 x 21.2 mm  
**Part No.:** [00B-4253-P0-AX](#)  
**Mobile Phase:** A: 0.1 % TFA in Water  
 B: 0.1 % TFA in Water / Acetonitrile (25:75)  
**Gradient:** Linear 93:7 (A/B) to 100% B over 5 minutes  
**Injection Volume:** 420 µL  
**Flow Rate:** 60 mL/min  
**Temperature:** Ambient  
**Detection:** UV @ 270 nm  
**Sample:** 1. Nadolol  
 2. Metoprolol  
 3. Propranolol


**Injection 240:** Axia Packed column with SecurityGuard PREP cartridge



**Time to change the SecurityGuard PREP cartridge**



**Injection 241:** Axia Packed column after removing SecurityGuard column protection system



**Original column performance maintained by using SecurityGuard PREP**





# Ordering Information



## SecurityGuard PREP

Semi-Prep	Preparative (HPLC or SFC)	
If your column ID (mm) is:		
9.0-16.0	18.0-29.0	30.0-49.0

Use cartridges (mm):		
10.0 x 10.0	15.0 x 21.2	15.0 x 30.0

Material	Description	pH Stability	10.0 x 10.0	15.0 x 21.2	15.0 x 30.0
<b>Cartridges for General Purpose/Pharmaceutical</b>					
<b>C18</b>	(ODS, Octadecyl)	1.5 - 10	AJO-7221	AJO-7839	AJO-8301
<b>C12</b>	(Dodecyl)	1.5 - 10	AJO-7275	AJO-7842	AJO-8304
<b>C8</b>	(MOS, Octyl)	1.5 - 10	AJO-7222	AJO-7840	AJO-8302
<b>C5</b>	(Pentyl)	1.5 - 10	AJO-7372	—	—
<b>Silica</b>	—	—	AJO-7223	AJO-7229	AJO-8312
<b>HILIC</b>	(HILIC)	1.5 - 8	AJO-8902	—	—
<b>NH<sub>2</sub></b>	(Amino, Aminopropyl)	1.5 - 11	AJO-7364	AJO-8162	AJO-8309
<b>CN</b>	(Cyano, Cyanopropyl)	2 - 7.5	AJO-7313	AJO-8220	AJO-8311
<b>Phenyl</b>	(Phenylhexyl)	1.5 - 10	AJO-7314	AJO-7841	AJO-8303
<b>PFP(2)</b>	(Pentafluorophenyl)	1.5 - 8	AJO-8376	AJO-8377	AJO-8378
<b>SCX</b>	(SA, Strong Cation Exchanger)	2.5 - 7.5	—	AJO-8595	AJO-8596
<b>SAX</b>	(SB, Strong Anion Exchanger)	2.5 - 7.5	AJO-7370	—	—
<b>RP-1</b>	(Reversed Phase - Polymer)	0 - 14	AJO-7368	AJO-8358	—
<b>Polar-RP</b>	(Ether-linked Phenyl)	1.5 - 7	AJO-7276	AJO-7845	AJO-8307
<b>Fusion-RP</b>	(C18 Polar Embedded)	1.5 - 10	AJO-7558	AJO-7844	AJO-8306
<b>AQ C18</b>	(Polar Endcapped C18)	1.5 - 7.5	AJO-7512	AJO-7843	AJO-8305
<b>Gemini<sup>®</sup> NX-C18</b>	(C18 TWIN-NX <sup>™</sup> Technology)	1 - 12	AJO-8369	AJO-8370	AJO-8371
<b>Gemini C18</b>	(C18 TWIN <sup>™</sup> Technology)	1 - 12	AJO-7598	AJO-7846	AJO-8308
<b>Gemini C6-Phenyl</b>	(C6-Phenyl TWIN Technology)	1 - 12	AJO-9156	AJO-9157	AJO-9158
<b>Luna<sup>®</sup> Omega Polar C18</b>	(Polar Function C18)	1.5 - 10	—	AJO-7603	AJO-7604
<b>Luna Omega PS C18</b>	(Mixed Mode C18)	1.5 - 10	—	AJO-7608	AJO-7609

**Cartridges for Core-Shell PREP Columns**

For core-shell media columns, such as Kinetex<sup>®</sup> and Aeris<sup>™</sup> (Phenomenex).

			/3 pk	/ea	/ea
<b>EVO C18</b>	(ODS, Octadecyl)	1 - 12	AJO-9306	AJO-9304	AJO-9305
<b>C18</b>	(ODS, Octadecyl)	1.5 - 8.5	AJO-9278	AJO-9145	AJO-9204
<b>C8</b>	(MOS, Octyl)	1.5 - 8.5	—	AJO-9205	AJO-9217
<b>PFP</b>	(Pentafluorophenyl)	1.5 - 9	—	AJO-9146	AJO-9299
<b>Phenyl-Hexyl</b>	(Phenylhexyl)	1.5 - 8.5	—	AJO-9147	AJO-9216
<b>Biphenyl</b>	(Biphenyl)	1.5 - 8.5	AJO-9280	AJO-9272	AJO-9273
<b>HILIC</b>	(HILIC)	2 - 7.5	—	AJO-9277	—
<b>C18-Peptide</b>	(ODS, Octadecyl)	1.5 - 9	AJO-9317	AJO-9318	AJO-9319
<b>F5</b>	(Pentafluorophenyl)	1.5 - 8.5	AJO-9323	AJO-9324	AJO-9325

**Cartridges for Protein and Polypeptide Reversed Phase**

For use with silica columns for separation of proteins and peptides, such as Jupiter<sup>®</sup> (Phenomenex); Vydac<sup>®</sup> 218TP, 214TP (Alltech Associates, Inc.); SynChropak<sup>®</sup> 300 C18, C4 (Eprogen, Inc.); Nucleosil<sup>®</sup> 300Å C18, C4 (Macherey-Nagel); Hypersil<sup>®</sup> 300Å (Thermo Hypersil-Keystone), and other widepore or 300Å brands.

			/3 pk	/ea	/ea
<b>Widepore C18</b>	(ODS, Octadecyl)	1.5 - 10	AJO-7224	AJO-7230	AJO-8313
<b>Widepore C5</b>	(Pentyl)	1.5 - 10	AJO-7371	—	—
<b>Widepore C4</b>	(Butyl)	1.5 - 10	AJO-7225	AJO-7231	AJO-8314

**Cartridges for Synthetic DNA / RNA Analysis**

For use with columns like Clarity<sup>®</sup> (Phenomenex)

			/3 pk	/ea	/ea
<b>Oligo-RP<sup>™</sup></b>	(C18 TWIN Technology)	1 - 12	AJO-8136	AJO-8210	—
<b>Oligo-WAX<sup>™</sup></b>	(WA, Weak Anion Exchanger)	1.5 - 11	AJO-8325	AJO-8639	—
<b>Oligo-XT</b>	(ODS, Octadecyl)	1 - 12	AJO-9516	AJO-9517	AJO-9518

**Cartridges for Silica GFC**

(Aqueous SEC) For use with silica GFC columns, such as Yarra<sup>™</sup> and BioSep<sup>™</sup> (Phenomenex); ZORBAX<sup>®</sup> GF-Series (Agilent Technologies); Bio-Sil<sup>®</sup> (Bio-Rad<sup>®</sup>).

				/ea	—
<b>GFC-2000</b>	—	2 - 7.5	—	AJO-8588	—
<b>GFC-3000</b>	—	2 - 7.5	—	AJO-8589	—
<b>GFC-4000</b>	—	2 - 7.5	—	AJO-8590	—

**Cartridges for Chiral**

For use with chiral columns, such as Lux<sup>®</sup> Cellulose-1, -2, -3, -4, i-Cellulose-5, and Amylose-1,-2 (Phenomenex); CHIRALCEL<sup>®</sup> OD-H<sup>®</sup>, CHIRALCEL<sup>®</sup> OJ-H<sup>®</sup>, and CHIRALPAK<sup>®</sup> AD<sup>®</sup>-H (DAICEL Corporation).

				/ea	/ea
<b>Lux i-Cellulose-5</b>	Cellulose tris(3,5-dichlorophenylcarbamate)	2 - 9	—	AJO-8634	AJO-8635
<b>Lux Cellulose-1</b>	Cellulose tris(3,5-dimethylphenylcarbamate)	2 - 9	AJO-8404	AJO-8405	AJO-8406
<b>Lux Cellulose-2</b>	Cellulose tris(3-chloro-4-methylphenylcarbamate)	2 - 9	AJO-8399	AJO-8400	AJO-8401
<b>Lux Cellulose-3</b>	Cellulose tris(4-methylbenzoate)	2 - 9	AJO-8623	AJO-8624	AJO-8625
<b>Lux Cellulose-4</b>	Cellulose tris(4-chloro-3-methylphenylcarbamate)	2 - 9	AJO-8628	AJO-8629	AJO-8630
<b>Lux Amylose-1</b>	Amylose tris(3,5-dimethylphenylcarbamate)	2 - 9	AJO-9344	AJO-9338	AJO-9339
<b>Lux Amylose-2</b>	Amylose tris(5-chloro-2-methylphenylcarbamate)	2 - 9	AJO-8472	AJO-8473	AJO-8474

**HPLC Guard Cartridge Holders\***

	/holder	/kit	/kit
<b>Reusable Holder</b>	AJO-9281	AJO-8223	AJO-8277

**SFC Guard Cartridge Holders\***

	/holder	/kit	/kit
<b>Reusable Holder</b>	AJO-9281	AJO-8617	AJO-8618



If SecurityGuard PREP cartridge protection system does not perform as well or better than your current guard cartridge system of similar phase and dimensions, return the product with the comparative data within 45 days for a FULL REFUND.

\*Includes column coupler

# Ordering Information

## Achiral Phases

### Aeris™

Phase	Length	ID	Part No.	Price
<b>5 µm</b>				
PEPTIDE XB-C18	150	21.2	<a href="#">00F-4632-P0-AX</a>	
	250	21.2	<a href="#">00G-4632-P0-AX</a>	

### Kinetex®

Phase	Length	ID	Part No.	Price
<b>5 µm</b>				
XB-C18	50	21.2	<a href="#">00B-4605-P0-AX</a>	
	50	30	<a href="#">00B-4605-U0-AX</a>	
	100	21.2	<a href="#">00D-4605-P0-AX</a>	
	100	30	<a href="#">00D-4605-U0-AX</a>	
	150	21.2	<a href="#">00F-4605-P0-AX</a>	
	150	30	<a href="#">00F-4605-U0-AX</a>	
	250	21.2	<a href="#">00G-4605-P0-AX</a>	
	250	30	<a href="#">00G-4605-U0-AX</a>	
EVO C18	50	21.2	<a href="#">00B-4633-P0-AX</a>	
	50	30	<a href="#">00B-4633-U0-AX</a>	
	100	21.2	<a href="#">00D-4633-P0-AX</a>	
	100	30	<a href="#">00D-4633-U0-AX</a>	
	150	21.2	<a href="#">00F-4633-P0-AX</a>	
	150	30	<a href="#">00F-4633-U0-AX</a>	
	250	21.2	<a href="#">00G-4633-P0-AX</a>	
	250	30	<a href="#">00G-4633-U0-AX</a>	
Biphenyl	100	21.2	<a href="#">00D-4627-P0-AX</a>	
	100	50	<a href="#">00D-4627-V0-AX</a>	
	150	21.2	<a href="#">00F-4627-P0-AX</a>	
	150	30	<a href="#">00F-4627-U0-AX</a>	
	250	21.2	<a href="#">00G-4627-P0-AX</a>	
HILIC	100	21.2	<a href="#">00D-4606-P0-AX</a>	
	150	21.2	<a href="#">00F-4606-P0-AX</a>	
	250	21.2	<a href="#">00G-4606-P0-AX</a>	
C18	50	21.2	<a href="#">00B-4601-P0-AX</a>	
	50	30	<a href="#">00B-4601-U0-AX</a>	
	100	21.2	<a href="#">00D-4601-P0-AX</a>	
	100	30	<a href="#">00D-4601-U0-AX</a>	
	150	21.2	<a href="#">00F-4601-P0-AX</a>	
	150	30	<a href="#">00F-4601-U0-AX</a>	
	250	21.2	<a href="#">00G-4601-P0-AX</a>	
	250	30	<a href="#">00G-4601-U0-AX</a>	
C8	50	21.2	<a href="#">00B-4608-P0-AX</a>	
	50	30	<a href="#">00B-4608-U0-AX</a>	
	100	21.2	<a href="#">00D-4608-P0-AX</a>	
	100	30	<a href="#">00D-4608-U0-AX</a>	
	150	21.2	<a href="#">00F-4608-P0-AX</a>	
	150	30	<a href="#">00F-4608-U0-AX</a>	
	250	21.2	<a href="#">00G-4608-P0-AX</a>	
	250	30	<a href="#">00G-4608-U0-AX</a>	
Phenyl-Hexyl	50	21.2	<a href="#">00B-4603-P0-AX</a>	
	50	30	<a href="#">00B-4603-U0-AX</a>	
	100	21.2	<a href="#">00D-4603-P0-AX</a>	
	100	30	<a href="#">00D-4603-U0-AX</a>	
	150	21.2	<a href="#">00F-4603-P0-AX</a>	
	150	30	<a href="#">00F-4603-U0-AX</a>	
	250	21.2	<a href="#">00G-4603-P0-AX</a>	
	250	30	<a href="#">00G-4603-U0-AX</a>	
F5	50	30	<a href="#">00B-4724-U0-AX</a>	
	100	30	<a href="#">00D-4724-U0-AX</a>	
	150	21.2	<a href="#">00F-4724-P0-AX</a>	
	150	30	<a href="#">00F-4724-U0-AX</a>	
	250	21.2	<a href="#">00G-4724-P0-AX</a>	
	150	30	<a href="#">00F-4603-U0-AX</a>	

**NEW**

### Gemini®

Phase	Length	ID	Part No.	Price
<b>5 µm</b>				
NX-C18	50	21.2	<a href="#">00B-4454-P0-AX</a>	
	50	30	<a href="#">00B-4454-U0-AX</a>	
	75	30	<a href="#">00C-4454-U0-AX</a>	
	100	21.2	<a href="#">00D-4454-P0-AX</a>	
	100	30	<a href="#">00D-4454-U0-AX</a>	
	150	21.2	<a href="#">00F-4454-P0-AX</a>	
	150	30	<a href="#">00F-4454-U0-AX</a>	
	250	21.2	<a href="#">00G-4454-P0-AX</a>	
	250	30	<a href="#">00G-4454-U0-AX</a>	
C18	50	21.2	<a href="#">00B-4435-P0-AX</a>	
	50	30	<a href="#">00B-4435-U0-AX</a>	
	100	21.2	<a href="#">00D-4435-P0-AX</a>	
	100	30	<a href="#">00D-4435-U0-AX</a>	
	150	21.2	<a href="#">00F-4435-P0-AX</a>	
	150	30	<a href="#">00F-4435-U0-AX</a>	
	250	21.2	<a href="#">00G-4435-P0-AX</a>	
	250	30	<a href="#">00G-4435-U0-AX</a>	
C6-Phenyl	100	21.2	<a href="#">00D-4444-P0-AX</a>	
	150	21.2	<a href="#">00F-4444-P0-AX</a>	
	250	21.2	<a href="#">00G-4444-P0-AX</a>	
<b>10 µm</b>				
NX-C18	50	21.2	<a href="#">00B-4455-P0-AX</a>	
	100	21.2	<a href="#">00D-4455-P0-AX</a>	
	100	30	<a href="#">00D-4455-U0-AX</a>	
	100	50	<a href="#">00D-4455-V0-AX</a>	
	150	21.2	<a href="#">00F-4455-P0-AX</a>	
	150	30	<a href="#">00F-4455-U0-AX</a>	
	150	50	<a href="#">00F-4455-V0-AX</a>	
	250	21.2	<a href="#">00G-4455-P0-AX</a>	
	250	30	<a href="#">00G-4455-U0-AX</a>	
	250	50	<a href="#">00G-4455-V0-AX</a>	
C18	100	21.2	<a href="#">00D-4436-P0-AX</a>	
	100	30	<a href="#">00D-4436-U0-AX</a>	
	150	21.2	<a href="#">00F-4436-P0-AX</a>	
	150	30	<a href="#">00F-4436-U0-AX</a>	
	150	50	<a href="#">00F-4436-V0-AX</a>	
	250	21.2	<a href="#">00G-4436-P0-AX</a>	
	250	30	<a href="#">00G-4436-U0-AX</a>	
	250	50	<a href="#">00G-4436-V0-AX</a>	

### Jupiter®

Phase	Length	ID	Part No.	Price
<b>4 µm</b>				
Proteo	250	30	<a href="#">00G-4396-U0-AX</a>	
<b>10 µm</b>				
Proteo	100	21.2	<a href="#">00D-4397-P0-AX</a>	
	250	21.2	<a href="#">00G-4397-P0-AX</a>	
	250	30	<a href="#">00G-4397-U0-AX</a>	
C18 300 Å	250	30	<a href="#">00G-4055-U0-AX</a>	
C4 300 Å	250	21.2	<a href="#">00G-4168-P0-AX</a>	

For additional sizes not displayed, please contact your Phenomenex technical consultant or local distributor.

guarantee

If Axia™ packed columns do not provide at least an equivalent separation as compared to a competing preparative column of the same particle size, same phase, and dimensions, return the product with comparative data within 45 days for a FULL REFUND. Only applies to 21.2 mm ID columns.



## Luna®

Phase	Length	ID	Part No.	Price
<b>5 µm</b>				
<b>C18(2)</b>	50	21.2	<a href="#">00B-4252-P0-AX</a>	
	50	30	<a href="#">00B-4252-U0-AX</a>	
	75	21.2	<a href="#">00C-4252-P0-AX</a>	
	75	30	<a href="#">00C-4252-U0-AX</a>	
	100	21.2	<a href="#">00D-4252-P0-AX</a>	
	100	30	<a href="#">00D-4252-U0-AX</a>	
	150	21.2	<a href="#">00F-4252-P0-AX</a>	
	150	30	<a href="#">00F-4252-U0-AX</a>	
	250	21.2	<a href="#">00G-4252-P0-AX</a>	
	250	30	<a href="#">00G-4252-U0-AX</a>	
<b>C8(2)</b>	75	30	<a href="#">00C-4249-U0-AX</a>	
	100	30	<a href="#">00D-4249-U0-AX</a>	
	150	21.2	<a href="#">00F-4249-P0-AX</a>	
	250	21.2	<a href="#">00G-4249-P0-AX</a>	
<b>CN</b>	250	21.2	<a href="#">00G-4255-P0-AX</a>	
<b>Phenyl-Hexyl</b>	150	21.2	<a href="#">00F-4257-P0-AX</a>	
<b>NH<sub>2</sub></b>	150	21.2	<a href="#">00F-4378-P0-AX</a>	
	250	21.2	<a href="#">00G-4378-P0-AX</a>	
<b>HILIC</b>	100	21.2	<a href="#">00D-4450-P0-AX</a>	
	150	21.2	<a href="#">00F-4450-P0-AX</a>	
	250	21.2	<a href="#">00G-4450-P0-AX</a>	
	250	30	<a href="#">00G-4450-U0-AX</a>	
<b>PFP(2)</b>	100	21.2	<a href="#">00D-4448-P0-AX</a>	
	100	30	<a href="#">00D-4448-U0-AX</a>	
	150	21.2	<a href="#">00F-4448-P0-AX</a>	
	250	21.2	<a href="#">00G-4448-P0-AX</a>	
	250	30	<a href="#">00G-4448-U0-AX</a>	
<b>Silica (2)</b>	100	21.2	<a href="#">00D-4274-P0-AX</a>	
	150	21.2	<a href="#">00F-4274-P0-AX</a>	
	250	21.2	<a href="#">00G-4274-P0-AX</a>	
	250	30	<a href="#">00G-4274-U0-AX</a>	
<b>10 µm</b>				
<b>C18(2)</b>	50	21.2	<a href="#">00B-4253-P0-AX</a>	
	100	21.2	<a href="#">00D-4253-P0-AX</a>	
	150	21.2	<a href="#">00F-4253-P0-AX</a>	
	150	30	<a href="#">00F-4253-U0-AX</a>	
	250	21.2	<a href="#">00G-4253-P0-AX</a>	
	250	30	<a href="#">00G-4253-U0-AX</a>	
	250	50	<a href="#">00G-4253-V0-AX</a>	
<b>C8(2)</b>	250	21.2	<a href="#">00G-4250-P0-AX</a>	
	250	50	<a href="#">00G-4250-V0-AX</a>	
<b>Silica (2)</b>	250	21.2	<a href="#">00G-4091-P0-AX</a>	
<b>15 µm</b>				
<b>C18(2)</b>	250	50	<a href="#">00G-4273-V0-AX</a>	
<b>C8(2)</b>	250	50	<a href="#">00G-4272-V0-AX</a>	

## Luna Omega

Phase	Length	ID	Part No.	Price
<b>5 µm</b>				
<b>Polar C18</b>	100	21.2	<a href="#">00D-4754-P0-AX</a>	
	100	30	<a href="#">00D-4754-U0-AX</a>	
	150	21.2	<a href="#">00F-4754-P0-AX</a>	
	150	30	<a href="#">00F-4754-U0-AX</a>	
	250	21.2	<a href="#">00G-4754-P0-AX</a>	
	250	30	<a href="#">00G-4754-U0-AX</a>	
	250	50	<a href="#">00G-4754-V0-AX</a>	
<b>PS C18</b>	50	21.2	<a href="#">00B-4753-P0-AX</a>	
	50	30	<a href="#">00B-4753-U0-AX</a>	
	100	21.2	<a href="#">00D-4753-P0-AX</a>	
	100	30	<a href="#">00D-4753-U0-AX</a>	
	150	21.2	<a href="#">00F-4753-P0-AX</a>	
	150	30	<a href="#">00F-4753-U0-AX</a>	
	250	21.2	<a href="#">00G-4753-P0-AX</a>	
	250	30	<a href="#">00G-4753-U0-AX</a>	
	250	50	<a href="#">00G-4753-V0-AX</a>	

Tip:

Protect your Axia™ Prep Columns with SecurityGuard™. See page 33 for ordering information

## Synergi™

Phase	Length	ID	Part No.	Price
<b>4 µm</b>				
<b>Fusion-RP</b>	100	21.2	<a href="#">00D-4424-P0-AX</a>	
	150	21.2	<a href="#">00F-4424-P0-AX</a>	
	250	21.2	<a href="#">00G-4424-P0-AX</a>	
<b>Hydro-RP</b>	50	21.2	<a href="#">00B-4375-P0-AX</a>	
	150	21.2	<a href="#">00F-4375-P0-AX</a>	
	250	21.2	<a href="#">00G-4375-P0-AX</a>	
<b>Max-RP</b>	150	21.2	<a href="#">00F-4337-P0-AX</a>	
	250	21.2	<a href="#">00G-4337-P0-AX</a>	
<b>Polar-RP</b>	50	21.2	<a href="#">00B-4336-P0-AX</a>	
	100	21.2	<a href="#">00D-4336-P0-AX</a>	
	100	30	<a href="#">00D-4336-U0-AX</a>	
	150	21.2	<a href="#">00F-4336-P0-AX</a>	
	150	30	<a href="#">00F-4336-U0-AX</a>	
	250	21.2	<a href="#">00G-4336-P0-AX</a>	
<b>10 µm</b>				
<b>Fusion-RP</b>	150	21.2	<a href="#">00F-4425-P0-AX</a>	
	250	21.2	<a href="#">00G-4425-P0-AX</a>	
<b>Hydro-RP</b>	150	21.2	<a href="#">00F-4376-P0-AX</a>	
	250	21.2	<a href="#">00G-4376-P0-AX</a>	
<b>Polar-RP</b>	250	21.2	<a href="#">00G-4351-P0-AX</a>	

## Clarity®

Phase	Length	ID	Part No.	Price
<b>5 µm</b>				
<b>Oligo-RP™</b>	100	21.2	<a href="#">00D-4442-P0-AX</a>	
	100	30	<a href="#">00D-4442-U0-AX</a>	
	250	21.2	<a href="#">00G-4442-P0-AX</a>	
<b>Oligo-XT</b>	100	21.2	<a href="#">00D-4745-P0-AX</a>	
	150	21.2	<a href="#">00F-4745-P0-AX</a>	
	150	30	<a href="#">00F-4745-U0-AX</a>	
	250	21.2	<a href="#">00G-4745-P0-AX</a>	
<b>10 µm</b>				
<b>Oligo-RP</b>	150	21.2	<a href="#">00F-4445-P0-AX</a>	
	150	30	<a href="#">00F-4445-U0-AX</a>	
	250	21.2	<a href="#">00G-4445-P0-AX</a>	
<b>Oligo-WAX™</b>	250	21.2	<a href="#">00G-4451-P0-AX</a>	

## Chiral Phases

### Lux®

Phase	Length	ID	Part No.	Price
<b>5 µm</b>				
<b>Amylose-1</b>	150	21.2	<a href="#">00F-4732-P0-AX</a>	
	250	21.2	<a href="#">00G-4732-P0-AX</a>	
	250	30	<a href="#">00G-4732-U0-AX</a>	
	250	50	<a href="#">00G-4732-V0-AX</a>	
<b>Amylose-2</b>	150	21.2	<a href="#">00F-4472-P0-AX</a>	
	250	21.2	<a href="#">00G-4472-P0-AX</a>	
	250	30	<a href="#">00G-4472-U0-AX</a>	
<b>Cellulose-1</b>	150	21.2	<a href="#">00F-4459-P0-AX</a>	
	250	21.2	<a href="#">00G-4459-P0-AX</a>	
	250	30	<a href="#">00G-4459-U0-AX</a>	
	250	50	<a href="#">00G-4459-V0-AX</a>	
<b>Cellulose-2</b>	150	21.2	<a href="#">00F-4457-P0-AX</a>	
	250	21.2	<a href="#">00G-4457-P0-AX</a>	
	250	30	<a href="#">00G-4457-U0-AX</a>	
	250	50	<a href="#">00G-4457-V0-AX</a>	
<b>Cellulose-3</b>	150	21.2	<a href="#">00F-4493-P0-AX</a>	
	250	21.2	<a href="#">00G-4493-P0-AX</a>	
	250	30	<a href="#">00G-4493-U0-AX</a>	
	250	50	<a href="#">00G-4493-V0-AX</a>	
	250	50	<a href="#">00G-4491-V0-AX</a>	
<b>Cellulose-4</b>	150	21.2	<a href="#">00F-4491-P0-AX</a>	
	250	21.2	<a href="#">00G-4491-P0-AX</a>	
	250	30	<a href="#">00G-4491-U0-AX</a>	
	250	50	<a href="#">00G-4491-V0-AX</a>	
	250	50	<a href="#">00G-4491-V0-AX</a>	
<b>i-Cellulose-5</b>	150	21.2	<a href="#">00F-4756-P0-AX</a>	
	250	21.2	<a href="#">00G-4756-P0-AX</a>	
	250	30	<a href="#">00G-4756-U0-AX</a>	
	250	50	<a href="#">00G-4756-V0-AX</a>	
	250	50	<a href="#">00G-4756-V0-AX</a>	

# The Ultimate Pre-Packed Preparative HPLC/SFC Column



## Australia

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auinfo@phenomenex.com

## Austria

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

## Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

## Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

## China

t: +86 400-606-8099  
f: +86 (0)22 2532-1033  
phen@agela.com

## Denmark

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

## Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

## India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

## Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

## Italy

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

## Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

## The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

## Norway

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

## Spain

t: +34 91-413-8613  
f: +34 91-413-2290  
espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

## USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

## All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com



## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at [www.Phenomenex.com/TermsAndConditions](http://www.Phenomenex.com/TermsAndConditions)

### Trademarks

Clarity, Gemini, Jupiter, Kinetex, Luna, and Lux are registered trademarks, and Aeris, BioSep, Oligo-RP, Oligo-WAX, Security-Guard, Synergi, TWIN-NX, TWIN Technology, Two-In-One Technology, and Yarra are trademarks of Phenomenex. Angiomax is a registered trademark of The Medicines Company. Bio-Sil and Bio-RAAd are registered trademarks of Bio-Rad Laboratories. NUCLEOSIL is a registered trademark of Macherey-Nagel. Hypersil is a registered trademark of Thermo Hypersil-Keystone LLC. SynChropak is a registered trademark of Eichrom Technologies, Inc. UPLC, XBridge, Symmetry, XTerra, and Waters are registered trademarks and OBD and SymmetryShield are trademarks of Waters Corporation. ZORBAX is a registered trademark of Agilent Technologies, Inc. Vydac is a registered trademark of Alltech Associates, Inc. CHIRALPAK, CHIRALCEL, AD, AD-H, AD-RH, AY, AY-H, OD, OD-H, OD-RH, OJ, OJ-H, OJ-RH, and OZ-H are registered trademarks of DAICEL Corporation.

### Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Waters Corporation, Agilent Technologies, Inc., Alltech Associates, Inc. or DAICEL Corporation. Columns used for comparison studies were manufactured by and purchased from Waters Corp., Agilent Technologies, Inc., Alltech Associates, Inc. and DAICEL Corporation.

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

Gemini and Kinetex EVO are patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

© 2017 Phenomenex, Inc. All rights reserved.